

# Eppendorf Certificate

## Multipette® M4/Repeater® M4

Multipette/Repeater M4 is the enhanced successor model of the Multipette/Repeater plus. Multipette/Repeater M4 does not differ in specification for any of the properties listed below.

- Technical specifications (accuracy + precision)
- Volume range from 1 µL to 10 mL
- Recognition of Combitips® size and display of dispensing volume
- Possibility to dispense 20 different volumes with one size of Combitip
- Positive displacement principle for dispensing
- Contamination free dispensing of hazardous liquids
- Precise & accurate dispensing of viscous solutions or liquids with high vapor pressure
- The optimization of the Multipette/Repeater M4 refers to the properties below.
- New step counter function
- Improved Combitip ejector enables one-handed Combitip ejection
- Optimized hand rest for enhanced ergonomics
- New sleep function for extending battery life
- Housing material with better chemical resistance
- New embedded RFID data chip with product properties
- Works only with improved Combitips advanced® for better tip attachment and more reliable wall dispensing
- Enlarged display size for a better readability
- Optimized tip sensor to detect a larger number of different tips and a better resistance against mechanical and chemical influences

Hamburg, April 2014



Christiane Markau  
Global Product Manager  
Liquid Handling



Dr. Birgit Schreiber  
Vice President  
Quality Management & Regulatory Affairs

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · 22331 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

ISO 9001  
Certified

ISO  
13485  
Certified

ISO  
14001  
Certified